




DESCRIPTION OF BREAST PUMPS









MANUAL BREAST PUMPS (SINGLE USER)


Reason for Issuance	<ul style="list-style-type: none"> - Attempting to breastfeed or provide breast milk. - Temporary issues such as engorgement, oversupply, plugged duct or sore nipples. - For mothers who need to pump for Infrequent separation from their baby, such as part-time return to work, school or other reason. 	
Qualifications	<ul style="list-style-type: none"> - NPP, BP or BE 	
PUMP MANUFACTURER, MODEL AND DESCRIPTION		
Ameda, Hand with Flexishield	The Flexishield Areola Stimulator™ provides increased stimulation and encourages the let-down reflex.	
Medela, Harmony	The dual pumping pattern that mimic a baby's sucking for the let-down and expression phases of breastfeeding. The PersonalFit Flex™ Breast Shield is designed with a new 105° angle and oval shape to reduce pressure on the breast and optimize milk flow.	
Hygeia, EnHande One Hand	Designed to ensure a direct flow from breast shield to milk storage container for maximum hygiene. One-piece valve system and fewer parts for easy assembly and cleaning. Silicone petal comfort massage shield helps with let-down.	

PERSONAL (SINGLE USER) ELECTRIC BREAST PUMPS

Reason for Issuance	<ul style="list-style-type: none"> - Experiencing separation for any period of time - Returning to work/school - Not appropriate for problems associated with initiating breastfeeding, poor latch, or temporary breastfeeding issues, such as engorgement or sore nipples.
Qualifications	<ul style="list-style-type: none"> - Categorized as exclusively breastfeeding (BE) and infant not receiving formula from WIC. - Established breastfeeding for at least one (1) month. - Client should receive only 1 personal use pump per pregnancy.

PUMP MANUFACTURER, MODEL AND DESCRIPTION

<p>Hygeia, EnJoye/EnRiche</p>	<p>The EnJoye is a cordless pump. The EnJoye and EnRiche mimic the infant’s unique suckling patterns with customizable speed and pressure controls. DISCONTINUED</p>	
<p>Hygeia Pro</p>	<p>The Hygeia Pro is a cordless pumps with a rechargeable battery, small and lightweight, simple to use, expression mode has 12 different settings, 13 different flange sizes available, ability to track pump sessions via app.</p>	
<p>Ameda, Finesse</p>	<p>Comfortflow™ technology, separate speed and suction dials with 32 fully customizable options for multiphase pumping. DISCONTINUED</p>	
<p>Medela, Personal Double</p>	<p>Single-user pump includes a double pump kit. Allows client to control suction strength and speed. DISCONTINUED</p>	
<p>Medela, Pump in Style</p>	<p>Single-user pump with a closed system, PersonalFit Flex breastshields. Includes battery pack</p>	
<p>Ardo Calypso Essential</p>	<p>Single-user pump with a closed system. Can operate with AA batteries. “Quietest pump on the market”. Shabbat Mode. DISCONTINUED</p>	
<p>Ardo Alyssa</p>	<p>Pump comes with a kit. Power pumping function that mimics cluster feeding. Allows saving of pump history with Memory Plus button.</p>	
<p>Mya Joy</p>	<p>Single-user pump. Travel size, compact, lightweight. Battery operated. Clips onto waistband.</p>	





Spectra S2 Plus	Personal pump weighing only 3lbs, the S2 Plus is ultra-quiet for optimum privacy. Additionally, the S2 Plus is hospital strength (270 mmHg), ensuring productive single or double pumping sessions. Like all Spectra® pumps, the S2 Plus is a true closed system, creating a hygienic barrier between outside contaminants and mother's milk.	
------------------------	---	--




MULTI-USER ELECTRIC BREAST PUMPS

Reason for Issuance	<ul style="list-style-type: none"> - Difficulty initiating and/or maintaining an adequate milk supply - Extended mother/infant separation - Maternal/infant illness - Desire to re-establish breastfeeding. - Other special circumstances determined by the Breastfeeding Coordinator.
Qualifications	<ul style="list-style-type: none"> - Categorized as NPP, BP or BE - Multi-user breast pumps should be issued for no longer than 1 month at a time. Situations requiring longer use must be re-evaluated at monthly intervals.

Each client will need their own attachment kit.

PUMP MANUFACTURER, MODEL AND DESCRIPTION

Ameda, Elite	Closed system eliminates the needs to clean and dry tubing. Able to control and adjust suction strength and cycling speed to mimic infant's suckling pattern.	
Ameda, Pearl	Two modes of pumping with independent speed and suction modes in expression with 72 settings. Lightweight at 3 lbs. Has a rechargeable battery. For use with Ameda Hygienikit.	
Hygeia, Endeare	Closed system eliminates the need to clean and dry tubing. Offers Customized Audio Recording Experience (CARE) for clients to record their infant's voice, cry or other sounds. Able to control and adjust suction strength and cycling speed to mimic infant's suckling pattern. DISCONTINUED	
Hygeia, Enjoye Cordless	multi-user pump and used with Universal Kit	

<p>Medela, Lactina</p>	<p>Allows for client to control suction strength. DISCONTINUED</p>	
<p>Medela, Symphony</p>	<p>Includes 2-Phase Expression that mimic a baby's sucking for the let-down and expression phases of breastfeeding.</p>	
<p>Calypso Pro</p>	<p>Multi-user version of the Calypso Essentials. AA Battery function</p>	

Developed by Oakland County Health Department. Revised October 2021 by Michigan Department of Health and Human Services WIC Division